



SIDE RAIL OPEN **ONLY** WHEN TRANSFERING BABY TO & FROM PARENT'S BED – CASTORS **MUST** BE LOCKED

CASTORS

KEEP ALL **FOUR** (4) CASTORS LOCKED AT ALL TIMES except when relocating





"up"



LOCKED "down"

Unlock ONLY when relocating



CLEANING

- Remove mattress
- Clean mattress, frame & shelf using manual scrubbing / wiping method
- Air dry in well ventilated area before replacing mattress

MAINTENANCE

- Wipe clean after each patient use
- Check for visible dirt / contamination & surface damage every day

Distributed in UK by:

